

Characterizing Glycan-Induced Structural Changes and Binding Interfaces for IgG1-C1q Complex using Hydroxyl Radical Protein Footprinting

Emily E Chea¹, Zhi Cheng¹, Jiana Duan¹, and Scot R Weinberger¹

1. GenNext Technologies, Inc., Half Moon Bay, CA, United States

HRPF Introduction

The Higher Order Structure (HOS) of a protein plays a critical role in a drug's stability, safety and biological function. Incorrect HOS or incorrect protein interactions are linked to adverse drug reactions which can result in further sickness or death. Advance techniques are required to robustly address the HOS of proteins. One such method is Hydroxyl Radical Protein Footprinting (HRPF). The Fox[®] Protein Footprinting System is a novel Hydroxyl Radical Protein Footprinting (HRPF) method that uses a proprietary flash oxidation lamp to generate hydroxyl radicals ($\cdot\text{OH}$) that irreversibly modify solvent exposed amino acid side chains. As solvent accessibility changes, the $\cdot\text{OH}$ modification concordantly changes.




Figure 1: Schematic of an HRPF method, fast photochemical oxidation of proteins (FPOP). With FPOP, protein is mixed with hydrogen peroxide and flown passed a pulsing light source which photolyzes the hydrogen peroxide into two $\cdot\text{OH}$ and modifies solvent exposed amino acids. Following labeling, the sample is deposited into a quench solution of catalase and methionine.

C1Q & IgG1 Introduction

The complement system is an essential part of our immune system that enhances the ability to clear pathogens and damaged cells, promotes inflammation, and attacks the pathogen's cell membrane. It consists of a series of proteins that work together in a cascade-like manner. The complement system can be activated via three pathways: the alternative pathway, the lectin pathway, and the classical pathway. The Classical pathway is initiated when the globular head of C1q which is apart of the C1 complex binds to antibodies that are attached to a pathogen. Once C1q's globular head binds to the antibody, it'll undergo a conformation change which will activate the C1 complex causing the continued cascade of enzymatic reactions ultimately leading to the lysis of the pathogen or infected cell.

C1q can bind to IgG and IgM. IgG is the most common type of antibody in blood circulation and is involved in response to most types of infections. There are four subtypes of IgG with IgG1 being the most common and it also has the highest binding affinity for C1q.

IgG has a single N-linked glycosylation site at asparagine 297 on each of the CH2 domains in the Fc region of the heavy chain. This glycosylation site is highly conserved in all four IgG subclasses. This glycan makes roughly 70 noncovalent interactions which stabilize the glycan and the Fc domain as well as how the Fab domain binds the antigen.

Previous studies have found IgG hexamerization is very important before C1q binds to the antibody. Further studies discovered the glycan on IgG is important for hexamer assembly which then results in C1q binding. However, there is still limited information in the full antibody characterization in native solution state conditions on exactly how the full antibody changes in structure following the removal of the glycan as well as how glycosylation effects C1q binding location and structure change.




Figure 2: Removing IgG1's glycan. Removed the glycan in native conditions of PBS and 37°C using PNGase F. We incubated at multiple time points, stopped the reaction by adding a denaturing buffer, digested with trypsin and then identified the peptide with the glycan. We searched for the peptides unmodified mass, deamidated mass, and the mass with the addition of the three most common glycans observed on this IgG1 (Adalimumab). After 3 hours we observed the majority of the glycans were removed with minor decrease with additional times.

Fox[®] Protein Footprinting System: HRPf labeling Method




Figure 3: Next we performed HRPf using (A) The Fox[®] System. All samples were incubated for six hours at 37 °C half with PNGase and half without. A single sample is 12 µL of IgG at 1.25 µM and about 3.5 minutes of labeling time. All samples and controls were acquired in triplicates. (B) The dose response plot for IgG1 with the glycan (green control and blue sample). (C) The dose response plot for IgG1 without the glycan. The dose response plots monitor the effective •OH concentration. This ensures we get the same change in absorbance/ •OH concentration between the conditions so any change in oxidation is a direct result from a change in interaction or structure. Following labeling the samples were digested with trypsin and were identified using Thermo's Vanquish Neo and Exploris 240. The data was processed using GenNext's FoxWare Software.

IgG1 ± Glycan HRPf Peptide Level Analysis

Glycosylation is known to stabilize the CH2 domain, however exactly how the structure changes has contrasting observations. Two x-ray crystallographic studies of de-glycosylated Fc showed the CH2 domain collapsing into a closed conformation while two other studies saw the opposite where the CH2 domain forms an open conformation. It is possible both take place, and the observed structures is highly dependent on the crystal packing environment. This points towards the importance of studying the antibody in native solution state conditions like HRPf. Not only does HRPf observe the more open flexible Fc domain but also provides information on the intact hinge region becoming protected and varying changes within the Fab region which is very important to keep in mind due to varying antigen binding following changes in glycosylation.




Figure 4: Peptide level HRPF Results. (A) Volcano plot comparing the $\text{Log}_2(\text{Fold Change})$ against $-\text{Log}(p\text{-value})$. Dotted line represents a p -value of 0.01. From this, we can easily identify multiple regions on IgG that change in average solvent accessibility following de-glycosylation. (B) Crystal structure of an IgG decorated to show regions with a significant change in solvent accessibility following the removal of the glycan. Regions in green were detected with no significant change in oxidation. Red peptides became more exposed while blue became more protected after glycan removal.

C1Q Peptide Level





Figure 5: Peptide Level HRPF results for C1Q \pm (IgG1 \pm glycan). (A) Histograms to show the total magnitude of oxidation for C1Q alone (dark blue), C1Q + IgG1-glycan (green), and C1Q + IgG1 + glycan (light blue). Peptides with a significant change in oxidation (p -value > 0.05) is marked with a red asterisks. (B) Volcano plot to show the fold change in oxidation and significance of change for C1Q with IgG1 + glycan (green) and C1Q with IgG1 - glycan (blue).

C1Q Residue Level





Figure 6: Residue Level HRPF results for C1Q \pm (IgG1 \pm glycan). (A) Histograms to show the total magnitude of oxidation for C1Q alone (dark blue), C1Q + IgG1-glycan (green), and C1Q + IgG1 + glycan (light blue). Peptides with a significant change in oxidation (p -value > 0.05) is marked with a red asterisks. (B) Volcano plot to show the fold change in oxidation and significance of change for C1Q with IgG1 + glycan (blue) and C1Q with IgG1 - glycan (green).

gC1Q Coverage Map

(A)	104	GS	GN	IK	DP	PP	AF	SA	RR	NP	PM	GG	NV	VI	FD	TV	IT	NQ	EE	PY	QN	HS	GR	FV	CT	VP	IG	YY	FT	F
	164	QV	LS	QW	EI	CL	SI	VSS	SR	GQ	VR	RS	LG	FD	TT	NK	GL	FQ	VV	SG	GM	VL	LQ	LQ	QD	QV	WV	EK	DP	K
	224	GH	IY	QG	SE	AD	SV	FS	GF	LI	FP	SA																		
(B)	118	AT	QK	IA	FS	AT	RT	IN	VP	LR	RD	QT	IR	FD	HV	IT	NM	NN	YE	PR	SG	KF	TC	KV	PG	LY	FT	YH	ASS	R
	178	GN	LC	VN	LR	GR	ER	AA	KV	VT	FC	DY	AY	NT	FQ	VT	TG	GM	VL	KL	LE	QE	EN	VF	LQ	AT	DK			
(C)	116	QK	FQ	SV	FT	LR	RQ	TH	DP	PA	NS	LR	FN	AA	TN	GM	GD	YD	TS	TG	KFT	CK	VP	GL	LY	FV	YH	AS	HT	
	176	AN	LC	VLL	LY	RS	GV	KV	VT	FC	GH	TS	KT	NQ	VNS	GG	VLL													




Figure 7: gC1Q coverage map for (A) subunit A (B) subunit B and (C) subunit C. Green font identifies peptide that were detected without a significant change in oxidation following peptide level analysis, red an increase in oxidation, and blue a decrease in oxidation. Residue cell's colored green were identified for residue level analysis with no change in oxidation, orange an increase in oxidation, and blue a decrease in oxidation.

Crystal structures decorated to identifies regions with a significant change in oxidation in (D) C1Q + IgG1-glycan and (E) C1Q + IgG1+glycan. The peptide backbone that is colored green are peptides identified with no significant change in oxidation following peptide level analysis, red an increase in oxidation, and blue a decrease. The residues that showed a significant change in oxidation are decorated as dot structures (red a significant increase in oxidation and blue a significant decrease in oxidation).

Conclusions

- HRPf identified glycan induced structure change for IgG1. The Fc region becomes more open, compacting around the hinge region resulting in multiple changes in the FAB region. The change in the FAB region points towards the significance of the glycan for antigen binding.
- HRPf identifies a binding pocket of IgG1 with C1Q with improved resolution using residues level analysis.
- C1Q binding decreases with de-glycosylated IgG1 but is not inhibited completely showing some similar regions of protection and exposure that are the same as IgG1 with the glycan. However, IgG1 with the glycan shows a larger fold change in oxidation.

GNTpNC1Q0923

Discover the Benefits
of Protein Footprinting

Contact us for products and services to investigate biopharmaceutical structure, interactions, folding, aggregation, formulation, and delivery.